

Read Alignment

Dr. Jung Soh
657.000 SS 2017

Mapping/Aligning reads

- ▶ Possible only when we have a reference
 - Reference genome or *de novo* assembly
- ▶ Issues
 - RNA-Seq reads spanning across exon junction
 - Alternative splicing
 - Reads mapping to multiple places in genome
- ▶ Different from multiple sequence alignment
 - MSA: align a group of sequences (similar length)
 - Read alignment: align reads to a region of reference sequence(s)

Why align reads to reference

- ▶ Assemble a new genome with reference alignment
- ▶ Analyze genetic variation from reference
- ▶ Check sequencing correctness
 - Most single-organism reads should align well to reference
- ▶ Analyze taxonomy of metagenomic reads
 - Align to multiple references
- ▶ Study differential expression of transcripts
 - Possible with RNA-Seq

Read alignment tools

- ▶ Many mapping tools
 - Bowtie, Bowtie2, SOAP, BWA, SHRiMP, mrFAST, mrsFAST, ZOOM, SSAHA2, Mosaik
- ▶ Mapping result
 - Most common format: SAM (sequence alignment/map)
 - Binary version: BAM
 - Use Samtools to analyze SAM/BAM files

Samtools

- ▶ Command line tool to work with SAM/BAM files
- ▶ Good for text-based analysis
 - SAM format contains lots of (often coded) information
 - Difficult to work with Linux commands
- ▶ Collection of commands
 - Select one command
 - Possibly provide options
 - Supply input SAM/BAM

End-to-end vs. local alignment

- ▶ **End-to-end: align all bases of a read**

Read: GACTGCGATCTCGACTTCG
Reference: TCGACTGGGCGATCTCGACTTCGAAAC
Alignment:

Read: GACTG--CGATCTCGACTTCG
 | | | | | | | | | | | | | | | | | | | | |

Reference: TCGACTGGGCGATCTCGACTTCGAAAC

- ▶ **Local: some bases at ends can be unaligned (clipped)**

Read: ACGGTTGCGTTAAATCCGCCACG
Reference: TAACTTGCGTTAAATCCGCCTGG
Alignment:

Read: ACGGTTGCGTTAA--TCCGCCACG
 | | | | | | | | | | | | | | | |

Reference: TAACTTGCGTTAAATCCGCCTGG

Bowtie2


- ▶ A widely used read alignment tool
- ▶ Two steps
 - Build index (database) from sequence files
 - Files can represent a genome, chromosome, or your own set of sequences
 - Align reads to the index
- ▶ Bowtie
 - Good for reads shorter than 50 bp
 - Only ungapped, end-to-end alignments
 - Read length upper limit of around 1000 bp

Bowtie2 key options

- ▶ Alignment mode
 - End-to-end (default)
 - Local
- ▶ Reporting policy
 - Search for multiple alignments, report the best one (default)
 - Search for up to N multiple alignments, report each
 - Search for and report all alignments

Alignment summary

- ▶ Usually given as alignment rate
- ▶ Can be complicated for paired-end reads
 - Especially with inconsistent terms used by alignment tools



Alignment visualization

- ▶ Similar to assembly visualization
 - Reads aligned to reference sequences (not assembled into contigs)
 - Reference sequences often annotated
- ▶ Many tools
 - IGV, SeqMonk, Tablet, BamView, Samtools